# YONESTEP

## DEUTSCH

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### FOR CULTIVATION OF HUMAN EMBRYOS

#### INDICATIONS FOR USE

YONESTEP is intended for the in vitro culture of human embryos following fertilization until day 5 to 6 of development.

#### PRODUCT DESCRIPTION

YONESTEP is a medium for culturing human embryos. It is suitable for cultivation of all stages of embryo development, for use with assisted reproduction techniques (ART).

#### COMPOSITION

- |                       |               |               |
|-----------------------|---------------|---------------|
| • Bicarbonate         | • Lactate     | • Glucose     |
| • Ultra Pure Water    | • Gentamicin  | • Phenol Red  |
| • Physiological Salts | • Hyaluronate | • Amino Acids |
| • Human Serum Albumin |               |               |

#### PRODUCT SPECIFICATIONS

- Endotoxin (LAL): < 0,25 EU/mL
- Sterility: sterile (no growth)
- pH: 7,20 - 7,40 (after equilibration in CO<sub>2</sub> atmosphere)
- MEA 1-cell to expanded blastocyst in 120h: ≥ 80%
- A certificate of analysis and MSDS are available upon request

**STERILE A** Sterilized using aseptic processing techniques.

#### KONTRAINDIKATION

No known contraindication.

#### WARNING

- Only for intended use.
- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discoloured, turbid or a precipitate is observed.
- Do not use if seal of the container is damaged upon delivery.
- Do not re-sterilize.

#### PRECAUTIONS

- R<sub>X</sub> Only** Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician.
- VITROMED recommends handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5).
  - Do not freeze before use.
  - Leave product in its original packaging.

- Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed colour.
- Discard minimal amounts of excess medium remaining in the bottle.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- YONESTEP contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods. Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, the product carries an extremely remote risk for transmission of viral diseases. The theoretical risk for the transmission of CJD is also considered extremely remote. No cases of viral disease or CJD have ever been identified for albumin.
- YONESTEP contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).
- The product may only be disposed of in accordance with the applicable regulations for waste products.

#### INSTRUCTIONS FOR USE

YONESTEP should be used only for its intended purpose by qualified and trained persons. VITROMED recommends that each laboratory uses its own optimized and validated laboratory procedures.

- 1) Before use equilibrate YONESTEP up to 37 °C and an atmosphere of 6% - 7% CO<sub>2</sub> in the CO<sub>2</sub> incubator for at least 4 hours, but optimally overnight. VITROMED recommends incubating aliquots. The medium must have contact with the CO<sub>2</sub> atmosphere during incubation (**Check the CO<sub>2</sub> content in the incubator with a measuring instrument and the pH in the medium after incubation!**).
- 2) Mix YONESTEP thoroughly before use, avoid foaming.
- 3) Remove the desired volume of product using aseptic techniques.
- 4) The medium is compatible with all culture flasks commonly used in the IVF laboratory ensuring gas exchange and can be used with or without mineral oil overlay.
- 5) VITROMED recommends using a 20 µL droplet containing one embryo. The culturing time is up to 144 h, a time lapse assessment system is recommended for embryo evaluation but not mandatory.
- 6) If using oil overlay add an amount of oil sufficient to cover the whole medium surface.

#### STORAGE INSTRUCTIONS AND STABILITY

- Do not use after expiry date.
- Store YONESTEP between 2 °C to 8 °C.
- The product should be used for a period of up to 7 days after opening, provided that it is handled under aseptic conditions.
- Keep away from sunlight.

#### REFERENCE

- Gardner Reprod Fertili Dev 2008; 20: 9-18.
- Reed et al. Fertili Steril 2009; 92: 1783-6.

## FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to [qc@vitromed.com](mailto:qc@vitromed.com). Report any serious incidents that have occurred with the product immediately by email to [qc@vitromed.com](mailto:qc@vitromed.com) and the competent authority of the Member State in which the user and/or patient is established or treated.

## SYMBOL GLOSSARY

(ISO 15223-1:2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements)

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Technischer Support  
Technical support

**VITROMED GmbH**  
Raiffeisenstr. 15a  
40764 Langenfeld  
Germany

+ 49 2173 . 200 41 30  
[support@vitromed.com](mailto:support@vitromed.com)  
[www.vitromed.com](http://www.vitromed.com)

Stand der Informationen / Date of information: 26-Nov-2020