## ENGLISH

### INDICATIONS FOR USE

V-VITWARM is intended for use in the thawing of vitrified of oocytes (MI), pronuclear (PN) zygotes through day 3 cleavage stage embryos and blastocyst stage embryos.

### CONTRAINDICATION

No known contraindication.

### PRODUCT DESCRIPTION

VITWARM is a ready-to-use kit to thaw human oocytes and embryos (up to blastocyst stage) which are cryopreserved. Subsequently, the thawed oocytes and embryos can be used in Assisted Reproduction Technology (ART) processes. There are in total 4 vials of pink solutions in the package:

- 2 vials of Thawing Solution (TS), with red caps and labels.
- 1 vial of Diluent Solution (DS), with a yellow cap and label.
- 1 vial of Washing Solution (WS), with a white cap and label.

### COMPOSITION

- Trehalose
- Hydroxypropyl cellulose
- Gentamicin

### QUALITY CONTROL TESTS PERFORMED

- Endotoxin (LAL, USP<85>): <0.25 EU/mL
- Sterility (USP<71>): sterile (no growth)
- pH (USP<791>): 7.20 – 7.45
- Osmolality (USP <785>):  
TS: 1250-1600 mOsm/kg (determined in 1:4 dilution)  
DS: 730-930 mOsm/kg (determined in 1:2 dilution)  
WS: 240-300 mOsm/kg (determined undiluted)
- One-cell MEA: ≥80% developed to expanded blastocyst at 96 hours
- Product appearance: solutions should be pink to rose color without precipitates

A certificate of analysis and MSDS are available upon request or can be downloaded from our website ([www.vitromed.com](http://www.vitromed.com)).

**STERILE A** Sterilized using aseptic processing techniques

### WARNING

- Do not use after expiry date.
- Do not freeze before use.
- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if seal the container (sterile barrier) is damaged upon delivery.
- Do not use medium that shows any evidence of particulate matter, cloudiness, or has changed colour.
- Do not re-sterilize.
- VITWARM contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).
- Note: The long-term safety of vitrification technique and maximum storage in liquid nitrogen has not been determined and unknown.

### PRECAUTIONS

- Prescription only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only for intended use.
- Read the instructions for use before use.
- This product is intended for use by medical specialists that have received fertility treatment training.
- Equipment and materials used should be sterilized.

- If eyes or skin come in contact with Vitrification/Thawing media, flush eye/skin with water immediately.
- Upon delivery media must be kept in original unopened container for storage and refrigerated at 2-8 °C.
- Before use, check the specifications of the container and labelling (cap color, vial labelling color, solution name, and volume). When finding something unusual in any of the above items, please quit using the product and get in touch with the distributor.
- VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5).
- Leave product in its original packaging.
- Excessive medium taken out of the primary package, which is not used up, must be discarded.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- The product may only be disposed of in accordance with the applicable regulations for waste products.

### PREPARATION BEFORE USE

- The whole process should be performed at room temperature unless otherwise indicated.
- Reverse the vials to mix the solutions before use. Shake the vials down using quick, sharp, downward wrist motions after reversing the vials to avoid liquid sticking in the inner side of the caps.
- Incubate the TS vial (closed) and the Petri dish in an incubator at 37 °C for at least 1.5 h.
- DS and WS should be at room temperature when thawing procedure starts.
- Use a sterilized pipette with appropriate diameter to handle oocytes (recommended 120 – 150 µm) and blastocysts (recommended 180 – 300 µm).
- Keep the cryopreserved oocytes or embryos frozen in liquid nitrogen, before the thawing procedure starts.
- Maximum 3-4 oocytes (MI) or pronuclear (PN) zygotes can be processed using these procedures before needing to replace medium.
- Maximum 1-2 embryos can be processed using these procedures before needing to replace medium.

### INSTRUCTIONS

- ① Take a sterile multi-well plate, and add 200 µl DS to the 1st well.
- ② Take the Petri dish TS vial out of the incubator, and add 1 ml to the Petri dish.
- ③ Take out one carrier device with a frozen oocyte/embryo from the cooling device (e.g., liquid nitrogen tank). (Note: Users in USA should use a closed cryopreservation device that is cleared in the US.) Within 1 sec, put the terminal of the device with the oocyte/embryo into the Petri dish with TS and incubate in the TS for 1 min. Make sure that the covering cap of the carrier device has been removed before placing the device into the TS. Make sure that the oocyte/embryo is immersed in the TS.
- ④ When in the TS, the oocyte/embryo naturally separates from the carrier device and start to float.
- ⑤ Aspirate the oocyte/embryo into a sterile pipette. After the oocyte/embryo is in the pipette, in addition, aspirate enough extra TS into the pipette (about 3 mm in length below the oocyte/embryo).
- ⑥ Insert the tip of the pipette into the centre bottom of the 1st well of the multi-well plate the with DS (prepared in step 1). Expel TS slowly to make a TS layer at the bottom of the well. Carefully place the oocyte/embryo into the bottom of the TS layer. Wait for 3 min.

- ⑦ During the waiting period, fill the 2nd and the 3rd wells with 200 µl WS.
- ⑧ Aspirate the oocyte/embryo in a sterile pipette, and then slowly aspirate extra DS in the pipette (about 3 mm in length below the oocyte/embryo).
- ⑨ Insert the tip of the pipette into the centre bottom of the 2nd well of the multi-well plate the with WS (prepared in step 7). Expel DS slowly to make a DS layer at the bottom of the well. Carefully place the oocyte/embryo into the bottom of the DS layer. Wait for 5 min.
- ⑩ After 5 min, observe the shape of the oocyte/embryo, if the shrunken oocyte/embryo returns to the original volume, it can be proceeded further.
- ⑪ Aspirate the oocyte/embryo in the pipette and transfer it to the 3rd well with WS. Place the oocyte/embryo at the surface of the WS and let the oocyte/embryo naturally sink to the bottom.
- ⑫ Take the oocyte/embryo again to the pipette, and place the oocyte/embryo at a different location of the surface of the WS in the same well. Let the oocyte/embryo naturally sink to the bottom.
- ⑬ Transfer the oocyte/embryo to a culture dish with appropriate culture medium. Assess the survival and proceed with further steps.

### STORAGE INSTRUCTIONS AND STABILITY

- The product is stable until the expiration date indicated on the label of each vial. Do not use, after the expiration date.
- Solutions must be stored at 2 – 8 °C.
- Keep away from sunlight.
- The product can be used for up to 7 days after opening the container, if it is handled under aseptic conditions.

### FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to [qc@vitromed.com](mailto:qc@vitromed.com). Report any serious incidents that have occurred with the product immediately by email to [qc@vitromed.com](mailto:qc@vitromed.com) and the competent authority of the Member State in which the user and/or patient is established.

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VITROMED GmbH