HYADASE

DEUTSCH

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ENGLISH

FOR ENZYMATIC REMOVAL OF CUMULUS AND CORONA RADIATA CELLS FROM OOCYTES

INDICATIONS FOR USE

HYADASE is intended for enzymatic removal of cumulus and corona radiata cells from oocytes prior to intracytoplasmic sperm injection (ICSI) procedures.

CONTRAINDICATION

No known contraindications.

PRODUCT DESCRIPTION

HYADASE is a ready-to-use media for the denudation of oocytes prior to the intracytoplasmic sperm injection (ICSI). It contains Hyaluronidase in HEPES buffered Human Tubal Fluid (HTF) Medium with Human Serum Albumin (HSA).

COMPOSITION

- | | | |
|---------------|-----------------------|--------------------|
| • HEPES | • Glucose | • Lactate |
| • Gentamicin | • Human Serum Albumin | • Ultra-Pure Water |
| • Bicarbonate | • Physiological Salts | • Amino acids |
| • Phenol Red | • Hyaluronidase | |

PRODUCT SPECIFICATIONS

- Endotoxin (USP<85>): < 1.0 EU/mL
- Sterility (USP <71>): sterile (no growth)
- Sterilized using aseptic processing techniques
- pH (USP<791>): 7.20 – 7.40
- Osmolality (USP<785>): 257 – 273 mOsm/kg
- MEA One-cell system: ≥80% of embryos developed to expanded blastocyst at 96 hours after a 4-minute exposure to V-HYADASE
- Hyaluronidase 50 - 120 IU/mL
- A certificate of analysis a MSDS and SSPC are available upon request.

WARNING

- Do not use after expiry date.
- Do not freeze before use.
- Do not re-sterilize.
- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discolored, turbid or a precipitate is observed.
- Do not use if the container(-seal) is damaged upon delivery.
- Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed color.
- HYADASE contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

PRECAUTIONS

- Prescription only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only for intended use.
- VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5).
- Leave product in its original packaging, reclose properly and store upright.
- Discard minimal amounts of excess medium remaining in the bottle.
- Discard media that has been warmed and not used.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- The product may only be disposed of in accordance with the applicable regulations for waste products.
- HYADASE contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods. Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, the product carries an extremely remote risk for transmission of viral diseases. The theoretical risk for the transmission of CJD is also considered extremely remote. No cases of viral disease or CJD have ever been identified for albumin.
- HYADASE contains pharmaceutical grade hyaluronidase which is derived from bovine testis and is certified with a TSE risk evaluation Certificate of Suitability (CEP). The animals from which the hyaluronidase is derived, are determined "fit for human consumption" and originate from countries with "negligible BSE risk", as determined Resolution "Recognition of the Bovine Spongiform Encephalopathy Risk Status of Member Countries", adopted by the OIE. According to the WHO guidelines on Tissue Infectivity Distribution in Transmissible Spongiform Encephalopathies (2010), testes from bovine source are classified as Category IC material (i. e. "Tissues with no detected infectivity").

INSTRUCTIONS FOR USE

HYADASE should be used only for its intended purpose by qualified and trained persons.

- 1 Prepare a petri dish containing one or more drops of HYADASE and several drops of VHEPES PLUS (80 µL to 100 µL each) and cover with VOIL. Equilibrate the dish to 37 °C for at least one hour.
- 2 Transfer up to 6 cumulus oocyte complexes (COC) into the HYADASE droplets and first gently pipette one COC up and down using a denuding pipette (e.g., VDENUKET) with a large inner diameter to remove oocytes from the cumulus cells.
- 3 After a maximum of 20 to 30 seconds transfer the partially denuded oocytes into a VHEPES PLUS droplet and remove the remaining cumulus or corona radiata cells by further denudation with a fine denuding pipette (e.g., VDENUKET).
- 4 Wash the oocytes thoroughly in the remaining VHEPES PLUS droplets.
- 5 Transfer the denuded oocytes in a dish with culture medium equilibrated for at least 4h in 37°C in a CO₂ incubator.
- 6 Oocytes should be left in culture medium for a short period of time (e.g., one hour) and are now ready for the ICSI procedure.

STORAGE INSTRUCTIONS AND STABILITY

- Store between 2 °C to 8 °C.
- Do not use after expiry date.
- The product can be used for a period of up to 7 days after opening, provided that it is handled under aseptic conditions.
- Keep away from sunlight.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com. Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the responsible competent authority.

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Technischer Support
Technical support

VITROMED GmbH
Raiffeisenstr. 15a
40764 Langenfeld
Germany

+49 21 73 .2 00 41 30
support@vitromed.com
www.vitromed.com

Stand der
Informationen
Date of information

17-April-2024

IFU-V-HYADASE_V03