DENUPET

DEUTSCH

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DENUPET

ENGLISH

FOR THE DENUDATION AND PIPETTING OF OOCYTES AND EMBRYOS

INDICATIONS FOR USE / INTENDED PURPOSE
[®]DENUPET is used for the removal of cumulus cells from cumulus oocyte complexes (COC) prior to Intracytoplasmic Sperm Injection (ICSI) and *In Vitro* Fertilization (IVF) as well as for the handling of oocytes and embryos during assisted reproductive techniques (ART). [®]DENUPET is not intended for biopsy of cells from oocytes or embryos.

PRODUCT DESCRIPTION
The tip of [®]DENUPET is available in a range of internal diameter sizes. Sizes 125 - 175 µm are suitable for oocyte denudation, sizes 200 - 600 µm are for oocyte and embryo handling procedures.

- PRODUCT SPECIFICATIONS**
- Material: Polycarbonate, non-pyrogenic
 - Total volumetric capacity (depending on the size): 20 – 25 µL
 - Proximal dimension: OD 900 µm
 - Endotoxins: < 0.05 EU/device
 - Sterility: sterile (SAL 10⁻⁶)
 - MEA 1-Cell to expanded blastocyst in 96h: ≥80%
 - A certificate of analysis and MSDS are available upon request.

STERILE R Sterilized by gamma irradiation

CONTRAINDICATION
No known contraindication.

- WARNINGS**
- Only for intended purpose.
 - Do not use if package is damaged.
 - Not sterile after opening.
 - Do not re-sterilize.
 - Do not reuse.
 - Multiple use may lead to in cross-contamination.
 - Each laboratory should decide for itself which [®]DENUPET diameter to use for each specific procedure! Internal diameter of the tip must be larger than the diameter of the sample. If the appropriate tip is not selected, sample damage is possible during ART procedures.
 - [®]DENUPET should be used only with pipettors compatible with the proximal dimensions of [®]DENUPET (0.9 mm).

- PRECAUTIONS**
- VITROMED recommends handling under aseptic conditions.
 - The products may only be disposed of in accordance with the applicable regulations for waste products.

Rx Only Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE
[®]DENUPET should only be used for its intended purpose by staff trained with the procedure.

- General handling:
- ① Attach a [®]DENUPET carefully on the pipetting device.
 - ② Rinse the [®]DENUPET by the pipetting device.
 - ③ Before working with the oocytes, embryos or blastocysts, aspirate a little medium into the [®]DENUPET. Only then start with the removal of the COC or the handling of oocytes and embryos.
 - ④ To avoid any air bubbles, do not completely empty fluid from the [®]DENUPET when using it.
 - ⑤ Dispose of the [®]DENUPET after use by appropriate measures.

- Denudation:
- ① Transfer the oocytes into media with an enzyme constituent, hyaluronidase, which promotes the denuding process.
 - ② Aspirate/expel one oocyte repeatedly through pipette tips of incrementally smaller internal diameters to remove the outer layers of cumulus and corona cells.
 - ③ Always keep the media level visible in the [®]DENUPET.

- STORAGE INSTRUCTIONS AND STABILITY**
- Do not use after expiry date.
 - Store at 15 – 25 °C.
 - Keep away from sunlight and heat.
 - Keep dry.

FEEDBACK AND OBLIGATION TO INFORM
Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com. Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established or treated.

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VITROMED GmbH