GRAD

DEUTSCH

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ENGLISH

INTENDED USE

GRAD is intended for separation of motile sperm from ejaculates for use in assisted reproduction procedures.

CONTRAINDICATION

No known contraindication.

PRODUCT DESCRIPTION

GRAD is a sterile colloidal suspension with silicate particles, stabilized with covalently bound hydrophilic silanes. Available as a 100 % stock solution and as a 40% or 80% ready-to-use solution with HEPES buffered Human Tubular Fluid (HTF) medium. GRAD provides highly efficient separation of motile sperm from most sperm samples.

COMPOSITION

- HEPES
- Silane-coated silica
- Gentamicin
- Phenol red

QUALITY CONTROL TESTS PERFORMED

- Endotoxin (LAL, USP<85>)
- Sterility (USP<71>)
- pH (USP<791>)
- Osmolality (USP<785>)
- Sperm survival test $\geq 80\%$ survival rate after 4 h exposure of sperm to test medium

A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.vitromed.com).

STERILE A Sterilized using aseptic processing techniques

WARNING

- Do not use after expiry date.
- Do not freeze before use.
- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if seal the container (sterile barrier) is damaged upon delivery.
- Do not use medium that shows any evidence of particulate matter, cloudiness, or has changed color.
- Do not re-sterilize.
- GRAD contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

PRECAUTIONS

- Prescription only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only for intended use.

- VITROMED recommends only handling of media with aseptic techniques (e. g., LAF-bench ISO Class 5).
- Leave product in its original packaging. Excessive medium taken out of the primary package, which is not used up, must be discarded.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- The product may only be disposed of in accordance with the applicable regulations for waste products.

PREPARATION BEFORE USE

- The whole process should be performed at room temperature.
- Use a sterilized pipette for withdrawing the desired volume of media
- Keep the cryopreserved sperms frozen at low temperature, before the thawing procedure starts.

INSTRUCTIONS

GRAD should only be used for its intended purpose, by trained personnel. VITROMED recommends the preparation of a two-step gradient system (45 % - 90 % or 40 % - 80 %) starting from GRAD 100%. If desired, a multilayer gradient (e.g. 45 % - 70 % - 90 %) can also be created.

1. Shake GRAD 100 % well before use.
2. To prepare a 90 % gradient medium, mix 1 part by volume of SPERMWASH medium with 9 parts by volume of GRAD 100 %.
3. To prepare a 45 % gradient medium, mix 5.5 parts by volume of SPERMWASH medium with 4.5 parts by volume of GRAD 100 %.

Note: Alternatively, other HEPES buffered media suitable for processing semen can be used to prepare the gradient media.

4. Mix the gradient media well to obtain homogeneous suspensions.
5. For separation of motile sperm from a semen sample, follow the instructions "GRADIENT SYSTEMS FOR SPERM SEPARATION".

Note: Gradient media should be prepared under sanitary conditions (LAF-bench ISO Class 5). For best results, prepare the gradient media no more than 24 hours before use.

GRADIENT SYSTEMS FOR SPERM SEPARATION

GRAD should only be used for its intended purpose, by trained personnel.

1. Shake GRAD well before use.
2. Warm all components and samples to room temperature.
3. Pour 2.5 ml of GRAD 80 % into a sterile disposable centrifuge tube.
4. Carefully layer 2.5 ml of GRAD 40 % on top of the GRAD 80 % without mixing the layers. Take care to maintain a clear phase boundary between the two layers.
5. Carefully pipette up to 2.5 ml of liquefied sperm onto the GRAD 40% layer.
6. Centrifuge for 15 - 18 minutes at 350 - 400 x g.

Note: If it is possible that no pellet can be seen. If this is the case, proceed with a second centrifugation for 3 - 5 minutes.

7. Carefully aspirate the supernatant down to 0.3 ml.
8. Resuspend the pellet in 2 - 3 ml of SPERMWASH, mix well and transfer to a new centrifuge tube.
9. Centrifuge the suspension at 300 x g for 8 - 10 min.
10. Aspirate the supernatant and repeat steps 8 and 9.
11. Aspirate the supernatant and resuspend the pellet in the desired volume of SPERMWASH.
12. The sample is now ready for further use.

Note: Highly viscous samples can be centrifuged up to a maximum of 500 x g. Avoid exceeding the specified centrifugal force.

STORAGE INSTRUCTIONS AND STABILITY

- The product is stable until the expiration date indicated on the label of each vial. Do not use after the expiration date.
- Solutions must be stored at 2 - 8 °C.
- Keep away from sunlight.
- The product can be used safely for up to 7 days after opening the container, if it is handled under aseptic conditions.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com.

Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established.

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* Ohne Phenolrot/ Without phenol red

VITROMED
Germany

Technischer Support
Technical support

VITROMED GmbH

Raiffeisenstr. 15a
40764 Langenfeld
Germany

+49 21 73 2 00 41 30
support@vitromed.com
www.vitromed.com

Stand der
Informationen
Date of information

18-Oct-2022

IFU-V-GRAD_V01