VGRAD 100, VGRAD 80, VGRAD 40, VGRAD KIT

DEUTSCH

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ENGLISH

INDICATIONS FOR USE

VGRAD is intended for separation of motile sperm from ejaculates by the density gradient method for use in assisted reproduction procedures.

CONTRAINDICATION

No known contraindication.

PRODUCT DESCRIPTION

VGRAD is a sterile colloidal suspension with silicate particles, stabilized with covalently bound hydrophilic silanes. Available as a 100 % stock solution (VGRAD 100) and as a 40% or 80% ready-to-use solution (VGRAD 40 and VGRAD 80) with HEPES buffered Human Tubular Fluid (HTF) medium. VGRAD 100 is a stock solution for preparing a density gradient system for semen preparation. VGRAD 40 and VGRAD 80 are gradient systems for semen preparation.

COMPOSITION

- HEPES
- Silane-coated silica
- Gentamicin
- Phenol red

PRODUCT SPECIFICATIONS

- Endotoxin (USP<85>): < 0.5 EU/mL
- Sterility (USP <71>): sterile (no growth)
- pH (USP<791>): 7.20 - 7.60
- Density (USP<841>): 1.050 - 1.061 g/cm³ (VGRAD 40)
1.100 - 1.111 g/cm³ (VGRAD 80)
1.125 - 1.136 g/cm³ (VGRAD 100)
- Osmolality: 310 - 340 mOsm/kg (VGRAD 40)
310 - 340 mOsm/kg (VGRAD 80)
300 - 330 mOsm/kg (VGRAD 100)
- HSSA: ≥ 80% of control motility at 24 hours after 1 hour exposure to test medium
- A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.vitromed.com).

STERILE

Sterilized using aseptic processing techniques

WARNINGS

- Do not use after expiry date.
- Do not freeze before use.
- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if seal the container (sterile barrier) is damaged upon delivery.
- Do not use medium that shows any evidence of particulate matter, cloudiness, or has changed color.
- Do not re-sterilize.
- VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5).
- VGRAD contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

PRECAUTIONS

- Prescription only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only for intended use.
- Leave product in its original packaging. Excessive medium taken out of the primary package, which is not used up, must be discarded.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- The product may only be disposed of in accordance with the applicable regulations for waste products.

INSTRUCTIONS FOR USE

VGRAD should only be used for its intended purpose, by trained personnel. The whole process should be performed at room temperature.

Preparation of gradients using VGRAD 100

VITROMED recommends the preparation of a two-step gradient system (40 % - 80 %) starting from VGRAD 100.

- ① Warm all components and samples to room temperature before use.
- ② Invert VGRAD 100 several times to mix the solution well before use.
- ③ To prepare a 80 % gradient medium, mix 1 part by volume of VSPERM WASH medium with 4 parts by volume of VGRAD 100.
- ④ To prepare a 40 % gradient medium, mix 6 parts by volume of VSPERM WASH medium with 4 parts by volume of VGRAD 100.
- ⑤ Mix the gradient media well to obtain homogeneous suspensions.
- ⑥ For separation of motile sperm from a semen sample, follow the instructions „GRADIENT SYSTEMS FOR SPERM SEPARATION“.

Note: Gradient media should be prepared under sanitary conditions (LAF-bench ISO Class 5). For best results, prepare the gradient media directly before use.

Gradient systems for sperm separation: Example using ready-to-use VGRAD 40 and VGRAD 80 (other gradient systems are possible)

- ① Warm all components and samples to room temperature before use.
- ② Invert VGRAD 40 and VGRAD 80 several times to mix the solutions well before use.
- ③ Pour 2.5 ml of VGRAD 80 into a sterile disposable centrifuge tube.
- ④ Carefully layer 2.5 ml of VGRAD 40 on top of the VGRAD 80 without mixing the layers. Take care to maintain a clear phase boundary between the two layers.

- ⑤ Carefully pipette up to 2.5 ml of liquified fresh or frozen sperm onto the °GRAD 40 layer.
- ⑥ Centrifuge for 15 - 18 minutes at 350 - 400 x g.
Note: It is possible that no pellet can be seen. If this is the case, proceed with a second centrifugation for 3 - 5 minutes.
- ⑦ Carefully aspirate the supernatant down to 0.3 ml.
- ⑧ Resuspend the pellet in 2 - 3 ml of °SPERM WASH, mix well and transfer to a new centrifuge tube.
- ⑨ Centrifuge the suspension at 300 x g for 8 - 10 min.
- ⑩ Aspirate the supernatant and repeat steps 8 and 9.
- ⑪ Aspirate the supernatant and resuspend the pellet in the desired volume of °SPERM WASH.
- ⑫ The sample is now ready for further use.
Note: Highly viscous samples can be centrifuged up to a maximum of 500 x g. Avoid exceeding the specified centrifugal force.

STORAGE INSTRUCTIONS AND STABILITY

- Store °GRAD between 2 - 8 °C.
- Keep away from sunlight.
- The product can be used for up to 7 days after opening the container, if it is handled under aseptic conditions.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to support@vitromed.com.

Report any serious incidents that have occurred with the product immediately by email to support@vitromed.com and the competent authority of the Member State in which the user and/or patient is established.

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* Ohne Phenolrot / Without phenol red

VITROMED
Germany

Technischer Support
Technical support

VITROMED GmbH
Raiffeisenstr. 15a
40764 Langenfeld
Germany

+49-2173-2004130
support@vitromed.com
www.vitromed.com

Stand der
Informationen
Date of information

15-Jan-2025

IFU-V-GRAD_V02