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ENGLISH

FOR THE COLLECTION, HANDLING, MICROMANIPULATION AND TRANSFER OF OOCYTES AND EMBRYOS

INDICATIONS FOR USE

HEPES PLUS is intended for oocyte collection, intracytoplasmic sperm injection (ICSI) fertilization, and oocyte and embryo washing and handling procedures. HEPES PLUS is also intended for transfer of embryos to the uterine cavity.

CONTRAINDICATION

No known contraindications.

PRODUCT DESCRIPTION

HEPES PLUS is a ready-to-use HEPES buffered Human Tubal Fluid (HTF) Medium.

COMPOSITION

- | | | |
|---------------|--------------------|-----------------------|
| • HEPES HTF | • Lactate | • Glucose |
| • Gentamicin | • Ultra Pure Water | • Human Serum Albumin |
| • Bicarbonate | • Amino acids | • Physiological Salts |
| • Phenol Red | | |

PRODUCT SPECIFICATIONS

- Endotoxin (LAL) (USP<85>): < 0.25 EU/mL
- Sterility (USP <71> / EP 2.6.1): sterile (no growth)
- pH: 7.00 - 7.40
- MEA 1-cell to expanded blastocyst in 96 h after 120-minute exposure: $\geq 80\%$
- A certificate of analysis and MSDS are available upon request.

STERILE A

Sterilized using aseptic processing techniques

WARNING

- Do not use after expiry date.
- Do not re-sterilize after opening.
- VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5).

- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discolored, turbid or a precipitate is observed.
- Do not use if seal of the bottle, or the bottle (sterile barrier) is damaged upon delivery.
- Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed color.
- HEPES PLUS contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).
- HEPES PLUS is not intended for follicle flushing procedures.

PRECAUTIONS

- Prescription only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only for intended use.
- Do not freeze before use.
- Leave product in its original packaging.
- Discard minimal amounts of excess medium remaining in the bottle.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- HEPES PLUS contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods. Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, the product carries an extremely remote risk for transmission of viral diseases. The theoretical risk for the transmission of CJD is also considered extremely remote. No cases of viral disease or CJD have ever been identified for albumin.
- The product may only be disposed of in accordance with the applicable regulations for waste products.

INSTRUCTIONS FOR USE

HEPES PLUS should be used only for its intended purpose by qualified and trained persons. VITROMED recommends that each laboratory uses its own optimized and validated laboratory procedures.

Oocyte Collection

- 1 Warm HEPES PLUS to 37 °C prior to use. Incubation should not be carried out in a CO₂ controlled incubator. As an alternative, use a screw-cap tube and ensure the cap is securely tightened during the incubation.
- 2 Remove the desired volume of product, sufficient to flush the equipment and for the procedure, using aseptic techniques. VITROMED recommends daily fresh aliquots in sterile tubes with well closing lids.
- 3 Place HEPES PLUS in sterile collection tubes to receive retrieved oocytes.
- 4 Use your own optimized and validated laboratory procedures for oocyte retrieval and follicle aspiration.
- 5 After collection, the oocytes can be transferred to a suitable culture medium for further IVF procedures or used for micro-manipulation.

Microinjection (ICSI)

- 1 Warm HEPES PLUS to 37 °C in culture dishes or flasks prior to use. Incubation should not be carried out in a CO₂ controlled incubator. As an alternative, use a screw-cap tube and ensure the cap is securely tightened during the incubation.

- 2
- Before microinjection place collected oocytes to remove cumulus cells into the wells or droplets on a dish loaded with hyaluronidase. After denudation incubate the oocytes in culture medium until microinjection.
- 3
- For the microinjection, transfer the denuded oocytes into VHEPES PLUS microdroplets overlaid with oil.
- 4
- Fertilize the oocyte by ICSI with prepared sperm and transfer the oocyte to an appropriate culture medium (e. g. VONESTEP).

Embryo Transfer
VHEPES PLUS is suitable for embryo transfer. The embryo can be placed in VHEPES PLUS for the transfer process. This method ensures a stable pH during transfer.

In certain practices, clinicians may choose to transfer embryos directly in the CO₂ pre-equilibrated bicarbonate-buffered culture medium, maintaining the embryo in its familiar environment. This approach supports continuity of conditions.

STORAGE INSTRUCTIONS AND STABILITY
• Do not use after expiry date.
• Store between 2 °C to 8 °C.
• The product can be used for a period of up to 7 days after opening, provided that it is handled under aseptic conditions.

FEEDBACK AND OBLIGATION TO INFORM
Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com.

Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established or treated.

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