INDICATIONS FOR USE

YVP is intended for use as an aid in the immobilization and isolation of individual sperm cells prior to intracytoplasmic sperm injection (ICSI) procedures.

CONTRAINDICATION

No known contraindications.

PRODUCT DESCRIPTION

YVP is a ready-to-use pharma grade polyvinylpyrrolidone (PVP) solution containing HEPES buffered Human Tubal Fluid (HTF) medium and Human Serum Albumin (HSA).

COMPOSITION

- | | | |
|------------------------|--------------|-----------|
| • Polyvinylpyrrolidone | • Lactate | • Glucose |
| • Human Serum Albumin | • Gentamicin | |
| • Physiological Salts | • HEPES HTF | |

PRODUCT SPECIFICATIONS

- Endotoxin (USP<85>): < 0.5 EU/mL
- Sterility (USP <71>): sterile (no growth)
- Sterilized using aseptic processing techniques
- pH: 7.25 – 7.45 (USP<791>)
- Osmolality: 300-330 mOsm/kg (USP<785>)
- Human Sperm Survival Assay ≥ 80% of control motility at 24h after 30 min contact to YVP.
- MEA (1-Cell System): ≥80% of embryos developed to expanded blastocyst at 96 h after 30-minute contact to YVP
- A certificate of analysis and MSDS are available upon request.

STERILE A Sterilized using aseptic processing techniques

WARNING

- Do not use after expiry date.
- Do not re-sterilize.
- Do not freeze before use.
- VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5).
- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discolored, turbid or a precipitate is observed.
- Do not use if seal of the container is damaged upon delivery.
- Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed color.

PRECAUTIONS

- Prescription only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only for intended use.
- Leave product in its original packaging.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- This product contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods. Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, the product carries an extremely remote risk for transmission of viral diseases. The theoretical risk for the transmission of CJD is also considered extremely remote. No cases of viral disease or CJD have ever been identified for albumin.

- YVP contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).
- Discard minimal amounts of excess medium remaining in the bottle.
- Discard media that has been warmed and not used.
- The product may only be disposed of in accordance with the applicable regulations for waste products.

INSTRUCTIONS FOR USE

YVP should be used only for its intended purpose by qualified and trained persons.

- ① Place several droplets of HEPES buffered media (5 µL – 10 µL) to hold prepared oocytes for the ICSI procedure and a droplet of YVP (5 µL – 10 µL) in a petri dish, cover with non-CO₂ preequilibrated YOIL and warm up to 37° C for at least 60 minutes.
- ② Add 1 µL – 2 µL of washed sperm to the edge of the YVP droplet.
- ③ Select, immobilize, and aspirate sperm on the opposite side and use it to inject the oocytes during ICSI.

STORAGE INSTRUCTIONS AND STABILITY

- Store YVP between 2° C to 8° C.
- The product can be used for a period of up to 7 days after opening, provided that it is handled under aseptic conditions.
- Keep away from sunlight.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com.

Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established or treated.

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VITROMED GmbH