DEUTSCH

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ENGLISH

INDICATIONS FOR USE

YOIL is intended for use as an overlay for culture of gametes, zygotes, or embryos in Assisted Reproduction Technology (ART) and micro-manipulation procedure.

INTENDED PURPOSE

YOIL is used to cover the medium during culturing and/or micromanipulation of human gametes and embryos.

PRODUCT DESCRIPTION

YOIL is a ready-to-use pre-washed, light paraffin oil.

CLINICAL BENEFIT

The clinical benefit of YOIL is related to its use in ART procedures for patients and couples with infertility problems. The overlay with YOIL prevents the evaporation of water and protects the media against changes in osmolality and pH. Culture of embryos in small volumes of medium under oil helps to enhance rates of embryonic development. The results of ART procedures using YOIL are expected to be comparable to ART outcomes (i.e. pregnancy rates, delivery rates) in the annual peer-reviewed report published by the European IVF-monitoring Consortium [EIM] for the European Society of Human Reproduction and Embryology [ESHRE].

COMPOSITION

- Mineral (paraffin) oil

PRODUCT SPECIFICATIONS

- Endotoxins (LAL) [USP<85>]: < 0,1 EU/mL
- Sterility [USP <71> / EP 2.6.1]: sterile (no growth)
- Density: 0,83 - 0,86 g/mL
- Viscosity: < 30 cP at 30 °C
- Peroxide Value (POV): < 0,1 mEq/kg
- MEA, 1-cell to expanded blastocyst in 120h: ≥80%
- Use of pharmaceutical grade mineral (paraffin) oil
- A certificate of analysis and MSDS are available upon request.

STERILE A Sterilized using aseptic processing techniques

CONTRAINDICATION

No known contraindication.

WARNING

- Do not use if the package is damaged.
- Not sterile after opening.
- Do not use if there are signs of microbial contamination.

- Do not use any medium that shows particles, turbidity or discoloration.
- Do not use if seal of the container is damaged upon delivery.

PRECAUTIONS

Rx Only Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

- Only for the intended use.
- VITROMED recommends handling media only with aseptic techniques (e.g. in a laminar flow workbench, ISO class 5).
- Leave product in its original packaging.
- Do not freeze before use.
- Do not reesterilize.
- Partial quantities that have been removed and not used must not be returned to the bottle.
- Minimal amounts of excess media in the bottle must be discarded.
- The products may only be disposed of in accordance with the applicable regulations for waste products.
- YOIL may appear cloudy when removed from the refrigerator or at room temperature after extended periods of storage time. Such cloudiness usually disappears when the oil is equilibrated in the incubator for several hours. If it does not disappear discard the product.
- Any human, organic material is considered potentially infectious. All samples should be handled as if they could transmit HIV or hepatitis.
- Always wear protective clothing when handling samples.

INSTRUCTIONS FOR USE

YOIL should be used only for its intended purpose by qualified and trained persons. Each laboratory should consult its own optimized and validated laboratory procedures.

- YOIL is pre-washed, no further preparation is necessary.
- VITROMED recommends preparing dishes/flasks with oil the day before use.
- Incubate dishes/flasks at 37°C in a CO₂ incubator so that the oil warms up and becomes saturated with gas.
- VITROMED recommends incubation of aliquots under aseptic conditions. If flasks are incubated, be sure not to tighten the cap of the flask because this will prevent gas exchange.
- Use a sterile, cell-culture suited pipette to add culture medium to a vessel which ensures gas exchange with the CO₂ atmosphere.
- Add an amount of oil sufficient to cover the whole medium surface.
- YOIL is compatible with all VITROMEDIA products.

STORAGE INSTRUCTIONS AND STABILITY

- Do not use after expiry date.
- Store V OIL between 2 °C to 8 °C.
- Keep away from sunlight.
- * V OIL can be used safely for at least 28 days after opening the container, provided that it is handled under aseptic conditions.
- Stable after transport (max. 5 days) at higher temperatures (≤ 37°C).


FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com.

Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established.

SYMBOL GLOSSARY

(ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by manufacturer - Part 1: General requirements)

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V-OIL-P100
V-OIL-P500
V-OIL-G100
V-OIL-G500

VITROMED
Germany

Technischer Support
Technical support

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VITROMED GmbH