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**ENGLISH****INTENDED USE**

VS<sup>®</sup>PERM FREEZE is used for cryopreservation and storage of human sperm.

**CONTRAINDICATION**

No known contraindication.

**PRODUCT DESCRIPTION**

VS<sup>®</sup>PERM FREEZE is a ready-to-use, HEPES buffered cryopreservation medium for use in *in vitro* fertilization (IVF) and intracytoplasmic sperm injection (ICSI) and is used for cryopreservation and storage of human sperm.

**COMPOSITION**

- HEPES
- Lactate
- Sucrose
- Glycine
- Glycerol
- Phenol red
- Human serum albumin (4 g/L)
- Physiological salts
- Dextrose monohydrate
- Gentamicin

**QUALITY CONTROL TESTS PERFORMED**

- Endotoxin (LAL, USP<85>)
- Sterility (USP<71>)
- pH (USP<791>)
- Osmolality (USP<785>)
- Sperm survival test  $\geq 80\%$  survival rate after 4 h exposure of sperm to test medium

A certificate of analysis and MSDS are available upon request or can be downloaded from our website ([www.vitromed.com](http://www.vitromed.com)).

**STERILE A**

Sterilized using aseptic processing techniques

**Warning**

- Do not use after expiry date.
- Do not freeze before use.
- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if seal the container (sterile barrier) is damaged upon delivery.
- Do not use medium that shows any evidence of particulate matter, cloudiness, or has changed color.
- Do not re-sterilize.
- VS<sup>®</sup>PERM FREEZE contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

**Precautions**

- Prescription only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only for intended use.
- VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5).
- Leave product in its original packaging.
- Excessive medium taken out of the primary package, which is not used up, must be discarded.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- The product may only be disposed of in accordance with the applicable regulations for waste products.

**PREPARATION BEFORE USE**

- The whole process should be performed at room temperature.
- Solutions can be warmed up to room temperature using a water bath or an incubator.
- Do NOT use any kinds of oven or microwave for the warming-up.
- Use a sterilized pipette for withdrawing the desired volume of media
- Keep the cryopreserved sperms frozen at low temperature, before the thawing procedure starts.

**INSTRUCTIONS****FREEZING**

- ① Liquefy the ejaculate for 30 minutes at room temperature before freezing.
- ② Mix 1 ml of semen with 1 to 0.7 ml of VS<sup>®</sup>PERM FREEZE. To do this, add medium dropwise to the semen and mix gently after each addition. Caution: To avoid cold shock, ensure that VS<sup>®</sup>PERM FREEZE is pre-warmed to room temperature.
- ③ Incubate for 10 minutes at room temperature.
- ④ Fill the sample mixture into a cryotube or cryo-straw according to the manufacturer's instructions.
- ⑤ For straws make sure that air bubbles are in the correct position according to the manufacturer's instructions.
- ⑥ Recommendation: Cool down the tubes or straws in a container with a low amount of water at 4°C for 20 minutes.
- ⑦ Dry the tube or straw and place it vertically in e.g., a Dewar vessel for 20 minutes, just above the level of liquid nitrogen.
- ⑧ Completely immerse cryotube or cryo-straw in the liquid nitrogen and transfer to a storage tank.

**DEFROST**

- ① Remove the desired cryotube or cryo-straw from liquid nitrogen and thaw for 10 minutes at room temperature. Then dry the straw off, so that no condensation water gets into the sample.
- ② Open the cryotube or cryo-straw according to the manufacturer's instructions and place the thawed sperm suspension in a sample tube.
- ③ Dilute sperm suspension with VS<sup>®</sup>PERMWASH (3 ml per 0.5 ml sperm), mix thoroughly and centrifuge at 300 x g to 350 x g for 15 minutes.
- ④ Resuspend the pellet in VS<sup>®</sup>PERMWASH.
- ⑤ Process the sperm suspension further by density gradient centrifugation or the swim-up method using VS<sup>®</sup>GRAD and VS<sup>®</sup>PERMWASH, depending on the use.

**STORAGE INSTRUCTIONS AND STABILITY**

- The product is stable until the expiration date indicated on the label of each vial. Do not use after the expiration date.
- Solutions must be stored at 2 – 8 °C.
- Keep away from sunlight.
- The product can be used safely for up to 7 days after opening the container, if it is handled under aseptic conditions.

**FEEDBACK AND OBLIGATION TO INFORM**

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com.

Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established.

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**VITROMED**  
Germany

Technischer Support  
Technical support

**VITROMED GmbH**  
Raiffeisenstr. 15a  
40764 Langenfeld  
Germany

+49 21 73 .2 00 41 30  
support@vitromed.com  
www.vitromed.com

Stand der  
Informationen  
Date of information

29-Sep-2022

IFU-V-SPERMFREEZE\_V01