VSPERMTAC

ENGLISH

INDICATIONS FOR USE

VSPERMTAC is intended for the determination of total antioxidant capacity in seminal plasma samples. For Research Purpose Only.

PRODUCT DESCRIPTION

VSPERMTAC 4 bottles dried reagent

COMPOSITION

- ABTS

REQUIRED MATERIALS (NOT INCLUDED THE KIT)

- VNACL0.9 – for reconstitution of the VSPERMTAC reagent and dilution of the sperm sample.
- PHOTOMETER – photometer for measuring absorbance according to the VSPERMTAC assay procedure.

PRODUCT SPECIFICATIONS

- Spectrophotometer test: passed
- Product appearance: passed
- A certificate of analysis and MSDS are available upon request.

WARNING

- Do not use after expiry date.
- Do not use the product if any signs of (microbial) contamination like decolorization, turbidity or a precipitate/cloudiness is observed.
- Do not use if the container is damaged upon delivery.
- Before first use of the Photometer, verify that the device is running the latest firmware version. If a firmware update is available, install it according to the manufacturer's instructions before operating the device.

PRECAUTIONS

- Only for intended use.
- Always wear protective clothing
- The product may only be disposed of in accordance with the applicable regulations for waste products.

INSTRUCTIONS FOR USE

VSPERMTAC should be used only for its intended purpose by qualified and trained persons.

- Open one bottle of VSPERMTAC and dissolve the dried reagent in 10 mL VNACL0.9. Incubate the reagent for 3 minutes and mix.
- Turn on the photometer and choose Program "CYA" (CYA = Capacity Yield of Antioxidants). Wait until the program is ready.
- Fill 10 mL VNACL0.9 solution in the middle cuvette and press "ZERO". Empty the measuring chamber.
- Fill the reagent in the middle measuring cuvette (2) of the Photometer, close the lid and press "OK TEST".
- Record the measured value (Value 1).
- Prepare the sperm dilution using VNACL0.9 according to the recommended volumes provided in the table below.

	Spermvolume [µl]	Diluent/VNACL0.9 [mL]
1	40	1
2	100	2,5
3	200	5

Using larger volumes reduces sample variability and increases pipetting accuracy. Take care to mix well.

- Transfer 200 µL of the diluted sperm sample into the middle cuvette (2). Mix thoroughly by pipetting up and down at least 5–10 times using a pipette with ≥5 mL capacity. Close the lid.
- Incubate 3 minutes

- Press "OK TEST"

- Record the second value (Value 2) and calculate the Antioxidant capacity of the sample according to the formula below, larger decrease indicates higher antioxidative potential.

$$\text{Antioxidant capacity [\%]} = \frac{\text{Value 1}}{\text{Value 2}} \cdot 100$$

- Immediately rinse the cuvette five times with water, and allow it to air dry. The instrument is waterproof; do not use cleaning agents.

STORAGE INSTRUCTIONS AND STABILITY

- Store between 2°C to 8°C.
- Use the reconstituted reagent within 2 hours after initial reconstitution.
- Protect from sunlight.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to support@vitromed.com.

Report any serious incidents that have occurred with the product immediately by email to support@vitromed.com.

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VITROMED GmbH