














SYMBOL GLOSSARY

Symbols according to EN ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by manufacturer - Part 1: General requirements

Symbol (reference number)	Title of symbol
 (5.1.1)	Manufacturer
 (5.1.4)	Use-by date
 (5.1.5)	Batch code
 (5.1.6)	Catalog number
 (5.2.5)	Steam sterilized
 (5.2.6)	Do not re-sterilize
 (5.2.8)	Do not use if package is damaged
 (5.2.11)	Single sterile barrier system
 (5.2.14)	Single sterile barrier system with protective packaging outside
 (5.3.1)	Fragile, handle with care
 (5.3.2)	Keep away from sunlight
 (5.3.4)	Keep dry
 (5.4.2)	Do not reuse

Date of information: 15.01.2025