

The background of the entire page features a large, semi-transparent petri dish on the left and a multi-well plate on the right. The petri dish has faint markings and numbers like '5' and '3'. The multi-well plate has a 'VITROMED Germany' label and numbered circles '1' and '3'.

IVF LABWARE

A VITROMED BRAND


Labware Catalogue

Our IVF Labware brand offers
secure solutions during your
laboratory procedures

www.ivflabware.com



IVF LABWARE of VITROMED is a specially designed portfolio of high quality medical plastic ware and particularly **IVF** tested for IVD use.



Focus on Quality Control

All our IVF Labware products are:

- **Cleanroom** manufactured
- **Individually packed**
- **Particle & scratch free**
- **Non toxic**
- **Endotoxin** tested
- 96h 1-cell Mouse Embryo Assay (**MEA**) > 80% tested
- **HSSA** > 70% tested
- Polymers with **USP Class VI** standard
- **Bisphenol A** free
- **Sterile by radiation**

V-CENTERWELL

Center-well-dish

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Center well Petridish	V-CWD-60	60 x 15	300



individual soft blister packed

V-CULTUREDISH-35, 60, 90

Petridishes with lid in different sizes

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Petridish 35 mm	V-CTD-35	35 x 10	300
Petridish 60 mm	V-CTD-60	60 x 15	300
Petridish 90 mm	V-CTD-90	90 x 20	300



individual soft blister packed

V-ICSIDISH-50

50 mm petridish

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
ICSI Petridish	V-ICD-50	50 x 10	300



individual soft blister packed

V-4 WELLDISH

4-well-dish

Details

- 4 well dish with flat bottom, volume 1.8 mL
- Wells marked with numbers
- Clear virgin polystyrene with excellent optical properties
- 200 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
4-Welldish	V-4WD	65 x 65	200



individual soft blister packed

V-ROUNDTUBE-5

Round bottom tube 5 mL

Details

- Dual position snap-cap
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box,
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Round bottom tube 5 ml	V-RTU-5	75 x 12	500



individual soft blister packed

V-ROUNDTUBE-14

Round bottom tube 14 mL

Details

- Dual position snap-cap
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Round bottom tube 14 ml	V-RTU-14	95 x 17	500



individual soft blister packed

V-CONICALTUBE-15

15 mL tube with conical shape

Details

- Printed graduations with white label area
- The practical screw cap prevents the liquid from leaking
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Centrifuge tube	V-CTU-15	120 x 17	500



individual soft blister packed

V-CRYOVIAL-2

2 mL vials for cryopreservation with screw cap

Details

- Round bottom, screw cap with white silicon ring
- Bottom specially designed for use in cryogenic racks and trays, free standing
- Polypropylene tubes to withstand temperatures to -196°C
- Printed graduations with white label area
- 50 pcs/bag, 1000 pcs/box

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Cryovial 2 ml	V-CRV-2	47 x 12	1000



V-SAMPLE CONTAINER

100 mL lockable sample container

Details

- Blue screw cap
- Graduations up to 90 mL, total volume 100 mL
- Clear virgin polypropylene with excellent optical properties
- 300 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Sperm samplecup	V-SAC-100	65 x 60	300



individual soft blister packed

V-SEROLOGICAL PIPET-1, 5, 10

Serological pipettes of different sizes

Details

- Printed +/- graduations
- color code for size identification
- Clear virgin polystyrene with excellent optical properties
- 300 & 200 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Serological pipet 1 ml	V-SRO-1	280 x 5	300
Serological pipet 5 ml	V-SRO-5	345 x 8	200
Serological pipet 10 ml	V-SRO-10	340 x 10	200



individual soft blister packed

V-TRANSFER PIPET-3

3 mL pasteur (transfer) pipette

Details

- Small tip ensures consistent drop size
- total volume of 3 mL
- Low density polyethylene
- 500 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Pasteur pipet 3 ml	V-TRP-3	150 length	500



individual soft blister packed

V-TIP-200, 1000

High quality pipette tips in different sizes

Details

- color code for size identification: 200 µl - yellow , 1000 µl - blue
- 500 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
200 µl pipet tip	V-TIP-200	n.a.	500
1000 µl pipet tip	V-TIP-1000	n.a.	500



individual soft blister packed





IVF

LABWARE

A VITROMED BRAND



VITROMED IVF LABWARE
manufacturing facility

JENA
Göschwitzer Str. 22
07745 Jena
Germany
Phone: +49 3641 539 19 76



other VITROMED
manufacturing facility:
LANGENFELD, Germany

other VITROMED offices:
Greece, Spain



VITROMED
Germany

www.ivflabware.com
order@vitromed.com

Connect with **VITROMEDonline**



For additional product or technical information please contact us at support@vitromed.com