

The background of the entire page features a large, semi-transparent petri dish on the left and a Vitromed IVF Labware tray on the right. The tray is labeled 'VITROMED Germany' and has circular wells numbered 1, 2, and 3. The text 'IVF LABWARE' is prominently displayed in the center, with 'A VITROMED BRAND' underneath it.

IVF LABWARE

A VITROMED BRAND

Labware Catalogue

Our IVF Labware brand offers
secure solutions during your
laboratory procedures

www.vitromed.com





Focus on Quality Control

All our IVF Labware products are:

- **Cleanroom** manufactured
- **Individually packed**
- **Particle & scratch free**
- **Non toxic**
- **Endotoxin** tested
- 96h 1-cell Mouse Embryo Assay (**MEA**) > 80% tested
- **HSSA** > 70% tested
- Polymers with **USP Class VI** standard
- **Bisphenol A** free
- **Sterile by radiation**



V-CENTERWELL

Center-well-dish

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Center-well-dish	V-CWD-60	60 x 15	300



individual soft blister packed

V-CULTUREDISH-35, 60, 90

Petridishes with lid in different sizes

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Petridish 35 mm	V-CTD-35	35 x 10	300
Petridish 60 mm	V-CTD-60	60 x 15	300
Petridish 90 mm	V-CTD-90	90 x 20	300



individual soft blister packed

V-ICSIDISH-50

50 mm petridish

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Petridish	V-ICD-50	50 x 10	300



individual soft blister packed

V-4 WELLDISH

4-well-dish

Details

- 4 well dish with flat bottom, volume 1.8 mL
- Wells marked with numbers
- Clear virgin polystyrene with excellent optical properties
- 200 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
4-Well-dish	V-4WD	65 x 65	200



individual soft blister packed

V-6 WELLDISH

6-well-dish

Details

- 6 well dish with flat bottom, volume 0,5 mL
- Clear, untreated virgin polystyrene with excellent optical properties
- High 12 mm, diameter of each well is 15 mm
- 200 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
6-Well-dish	V-6WD	65 x 75 mm	200



individual soft blister packed

V-ROUNDTUBE-5

Round bottom tube 5 mL

Details

- Dual position snap-cap
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Round bottom tube 5 ml	V-RTU-5	75 x 12	500



individual soft blister packed

V-ROUNDTUBE-14

Round bottom tube 14 mL

Details

- Dual position snap-cap
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Round bottom tube 14 ml	V-RTU-14	95 x 17	500



individual soft blister packed

V-CONICALTUBE-15

15 mL tube with conical shape

Details

- Printed graduations with white label area
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Centrifuge tube	V-CTU-15	120 x 17	500



individual soft blister packed

V-CRYOVIAL-2

2 mL vials for cryopreservation with screw cap

Details

- Round bottom, screw cap with white silicon ring
- Bottom specially designed for use in cryogenic racks and trays, free standing
- Polypropylene tubes to withstand temperatures to -196°C
- Printed graduations with white label area
- 50 pcs/bag, 1000 pcs/box

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Cryovial 2 ml	V-CRV-2	47 x 12	1000



V-SAMPLE CONTAINER

100 mL lockable sample container

Details

- Blue screw cap
- Graduations up to 90 mL, total volume 100 mL
- Clear virgin polypropylene with excellent optical properties
- 300 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Sperm samplecup	V-SAC-100	65 x 53 (without lid)	300



individual soft blister packed

V-SEROLOGICAL PIPET-1, 5, 10

Serological pipets of different sizes

Details

- Printed +/- graduations
- Color code for size identification
- Clear virgin polystyrene with excellent optical properties
- 300 & 200 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Serological pipet 1 ml	V-SRO-1	280 x 5	300
Serological pipet 5 ml	V-SRO-5	345 x 8	200
Serological pipet 10 ml	V-SRO-10	340 x 10	200



individual soft blister packed

V-TRANSFER PIPET-3

3 mL pasteur (transfer) pipet

Details

- Small tip ensures consistent drop size
- total volume of 3 mL
- Low density polyethylene
- 500 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Pasteur pipet 3 ml	V-TRP-3	150 length	500



individual soft blister packed

V-TIP-200, 1000

High quality pipet tips in different sizes

Details

- color code for size identification: 200 µL- yellow , 1000 µL- blue
- 500 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
200 µl pipet tip	V-TIP-200	n.a.	500
1000 µl pipet tip	V-TIP-1000	n.a.	500



individual soft blister packed

V-HANDPET

Plastic pipet for liquid handling

Details

- Made of polycarbonate tubing
- Individual soft blister packed
- Radiation sterilized
- Compatible with luer lock syringes

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Plastic hand pipet	V-HPT-90	90	250



individual soft blister packed



IVF

LABWARE

A VITROMED BRAND



VITROMED IVF LABWARE
manufacturing facility

JENA
Göschwitzer Str. 22
07745 Jena
Germany
Phone: +49 3641 539 19 76



other VITROMED
manufacturing facility:
LANGENFELD, Germany

other VITROMED offices:
Greece, Spain



VITROMED
Germany

www.vitromed.com
order@vitromed.com

Connect with **VITROMEDonline**



For additional product or technical information please contact us at support@vitromed.com