

Product Catalogue

We stand for German precision, high quality products and compassionate customer service

www.vitromed.com



Who we are

VITROMED was founded in 2011 and is now a leading global company in manufacturing high-quality medical devices used during all stages in fertility treatment in IVF clinics. Our advanced medical device manufacturing facilities are located in Germany, Jena and Langenfeld, but our products are available worldwide. In 2020 and 2021 VITROMED established representative offices in Spain and Greece for marketing and other service operations.

Our Commitment to Quality & Service

VITROMED is an ISO 13485 certified company that ensures the manufacturing of top-notch products through a quality management system and with the usage of advanced technologies and the latest research. Most VITROMED products are MEA and LAL tested and undergo stringent, lot-wise quality control before product release.

Most of VITROMED products are CE marked and FDA registered. More information can be found on our website.

We strive to supersede expectations and take extreme care in meeting customer requirements. Over the years, we have been providing support with authentic products and services to numerous clinics around the globe. We ensure our customer satisfaction with high-quality customized products augmented with beautiful design and integrity.

See what our clients say about us at www.vitromed.com

Ask our Scientific Advisors

We have well-trained and highly professional technical and service teams, who are always ready to resolve any problem our customers face and provide information regarding our products.

Additionally, we have a team of experts, headed by our Scientific Advisors, the late Dr. Patrick Quinn and Dr. Ralf Böhm, who are leading several research and development projects.

VITROMEDIA —	4
MICROTOOLS	8
IVD & ACCESSORIES —	10
V-CARE	15



VITROMEDia



GRAD



V-GRAD is intended for separation of motile sperm from ejaculates by the density gradient method for use in assisted reproduction procedures.

Details

- Endotoxin tested (< 0.5 EU/mL)
- Human Sperm Survival Assay tested ≥ 80% Control Motility
- Colored with Phenol Red (V-GRAD 40 and 80)
- V-GRAD 100% stock solution
- Ready-to-use V-GRAD 40% and 80% diluted with HEPES buffered HTF
- Store at 2-8 °C
- 1 x 100 mL PETG bottle and glass bottle
- Kit: 2 x 12 mL V-GRAD 40 and 2 x 12 mL V-GRAD 80 glass bottles

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
Gradient 40	V-GRD43_100P	100	PETG	1
Gradient 80	V-GRD83_100P	100	PETG	1
Gradient 100	V-GRD00_100P	100	PETG	1
Gradient Kit	V-GRDK4_12	4x12	GLASS	1 Kit





 $V\hbox{-SPERM}\, FREEZE is a ready-to-use \, HEPES \, buffered \, medium \, and \, used \, for \, the \, freezing \, of \, human \, sperm.$

Details

- Endotoxin tested (< 0.25 EU/mL)
- Human Sperm Survival Assay tested ≥ 80% Control Motility
- Store at 2-8 °C
- 4 x 12 mL glass bottles

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
Sperm Freeze	V-SFM-K2_12G	12	GLASS	4 x 12





V-OIL is a ready to use light paraffin oil for use as an overlay for culture of gametes, zygotes, or embryos in Assisted Reproduction Technology (ART) and micromanipulation procedure.

Details

- Endotoxin tested (< 0.1 EU/mL)
- MEA 1-cell to expanded blastocyst in 120h: \geq 80%
- Blastocyst cell count on day 6 of mouse embryo development
- Peroxide value (POV) determination (<0.1 mEq/KG)
- Prewashed twice with high-purity water
- Available in 100 mL and 500 mL with PETG/PET and glass bottles

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
Paraffin Oil	V-OIL-P100	100	PETG	1
Paraffin Oi l	V-OIL-P500	500	PETG	1
Paraffin Oil	V-OIL-G100	100	GLASS	1
Paraffin Oil	V-OIL-G500	500	GLASS	1





FDA Registered

V-SPERM WASH is a HEPES buffered HTF medium with HSA, intended for washing of spermatozoa, the isolation of motile viable sperm by swim-up method and holding sperm for IUI procedures.

Details

- Endotoxin tested (< 0.25 EU/mL)
- Human Sperm Survival Assay tested ≥ 80% Control Motility
- Store at 2-8 °C
- 1 x 100 mL PETG/PET bottle

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
Sperm Wash	V-SWM-100	100	PET(G)	1





FDA Registered

V-PVP is intended for use as an aid in the immobilization and isolation of individual sperm cells prior to intracytoplasmic sperm injection (ICSI) procedures.

Details

- Endotoxin tested (< 0.5 EU/mL)
- Human Sperm Survival Assay tested ≥ 80% Control Motility
- MEA 1-Cell to expanded blastocyst in 96h: ≥ 80%
- Store at 2-8 °C
- 5 x 0.5 mL Microvials

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
Sperm Immobilization	V-PVP-7	0.5	PPCO	5 x 0.5 mL
Sperm Immobilization	V-PVP-10	0.5	PPCO	5 x 0.5 mL

HYADASE



FDA Registered

 $V-HYADASE\ is\ intended\ for\ enzymatic\ removal\ of\ cumulus\ and\ corona\ radiata\ cells\ from\ oocytes\ prior\ to\ intracytoplasmic\ sperm\ injection\ (ICSI)\ procedures.$

Details

- 50 -120 U /ml Hyaluronidase
- Endotoxin tested (< 1.0 EU/mL)
- MEA 1-Cell to expanded blastocyst in 96h: \geq 80%
- Easy and secure removal of cumulus and corona radiata cells

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
Hyaluronidase	V-HYA	1	PPCO	5 x 1





FDA Registered

 $V-HEPES\ PLUS\ is\ intended\ for\ oocyte\ collection, intracytoplasmic\ sperm\ injection\ (ICSI)$ fertilization, and oocyte\ and\ embryo\ washing\ and\ handling\ procedures. $V-HEPES\ PLUS\ is\ also\ intended\ for\ transfer\ of\ embryos\ to\ the\ uterine\ cavity.$

Details

- Endotoxin tested (< 0.25 EU/mL)
- MEA 1-Cell to expanded blastocyst in 96h: ≥ 80%
- Store at 2-8°C
- 1 x 100 mL PETG /PET bottle

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
HEPES PLUS	V-HTF-P1	100	PET(G)	1





FDA Registered

 $V-ONESTEP\ is\ intended\ for\ the\ in\ vitro\ culture\ of\ human\ embryos\ following\ fertlization\ until\ day\ 5\ to\ 6\ of\ development.$

Details

- Endotoxin tested (< 0.25 EU/mL)
- MEA 1-cell to hatching blastocyst in 120h: ≥80%
- Blastocyst cell count on day 6 of mouse embryo development
- Time-lapse morphokinectis MEA for 120 h
- Store at 2-8 °C
- 1 x 20 mL Glass, 1 x 20 mL screw cap PETG, 1x 10 ml Glass, 1 x 10 mL screw cap PETG and 1 x 50 mL screw cap PETG bottle

Ordering Information*

Ref Nr.	Volume (mL)	Packaging	Qty/Unit
V-OSM-20	20	GLASS	1
V-OSM-20P	20	PETG	1
V-OSM-10G	10	GLASS	1
V-OSM-10	10	PETG	1
V-OSM-60	50	PETG	1
	V-OSM-20 V-OSM-20P V-OSM-10G V-OSM-10	V-OSM-20 20 V-OSM-20P 20 V-OSM-10G 10 V-OSM-10 10	V-OSM-20 20 GLASS V-OSM-20P 20 PETG V-OSM-10G 10 GLASS V-OSM-10 10 PETG

^{*}check the product avability in your country

Info

Time-lapse technology is a perfect solution for observing the embryo development while culturing without removing the embryos from incubators. This method offers much more details on the timing of the cleavages and the embryo morphology. V-OIL and V-ONESTEP undergo time-lapse morphokinectis MEA for 120 h to provide the most sophisticated QC results.

Please check the product availability in your region while ordering.

MICROTOOLS



ICSIPET



V-ICSIPET HOLDING and INJECTION pipettes are manufactured from special cleaned borosilicate glass under cleanroom conditions.

Details

- Endotoxin tested (< 0.05 EU/device)
- MEA 1-Cell to expanded blastocyst in 96h: ≥ 80%
- Gamma sterilized
- Customized single packaging for easy and secure removal
- Color-coded according to purpose
- 1 x 10 pipettes individually, sterile packed

Ordering Information

Туре	Ref Nr.	Dimension (µm)	Angle (°)	Qty/Unit
INJECTION	V-INJ-S3-30	7 (OD)	30	10
INJECTION	V- I NJ-S3-35	7 (OD)	35	10
HOLDING	V-HOLD-30	120 (OD)	30	10
HOLDING	V-HOLD-35	120 (OD)	35	10

DENUPET



 $V\hbox{-}DENUPET\ pipettes\ are\ manufactured\ from\ high\ quality\ medical\ polymer\ under\ clean room\ conditions.$

Details

- Endotoxin tested (< 0.05 EU/device)
- MEA 1-Cell to expanded blastocyst in 96h: ≥ 80%
- Gamma sterilized
- Easy and secure removal
- Unique color code for each size V-DENUPET - 5 x 10 pipettes

Ordering Information

Туре	Ref Nr.	Dimension (μm)	Qty/Unit
DENUPET	V-DEN-	125, 135, 140, 150, 175 200, 275, 300, 600	50





Biopsy pipettes are manufactured from special cleaned borosilicate glass under cleanroom conditions.

Details

- Endotoxin tested (< 0.05 EU/device)
- MEA 1-Cell to expanded blastocyst in 96h: ≥ 80%
- Gamma sterilized
- Customized single packaging for easy and secure removal
- Color-coded according to purpose
- $\bullet \quad 1 \times 10$ pipettes individually, sterile packed

Ordering Information

Туре	Ref Nr.	Dimension (µm)	Angle (°)	Qty/Unit
Polar Body	V-PGDP-V1-30	23 (OD)	30	10
Polar Body	V-PGDP-V1-35	23 (OD)	35	10
Blastomere*	V-BBP-V1-30	46 (OD)	30	10
Blastomere	V-BBP-V1-35	46 (OD)	35	10
Trophectoderm*	V-TEBP-V1-30	32 (OD)	30	10
Trophectoderm	V-TEBP-V1-35	32 (OD)	35	10

For custom dimensions please contact our support team *FDA Registered only

Info

Please check the product availability in your region while ordering.

IVD & Accessories





V-DENUPET HANDLE is used in combination with denuding pipettes during assisted reproductive treatment for removement of cumulus cells from cumulus-oocyte complexes (COC), handling of oocytes, pronuclear stages, embryos and blastocysts.

Details

- V-DENUPET HANDLE allows secure volume adjustments till 3 μ L
- Precise control for manipulation of oocytes and embryos
- Smooth and responsive plunger making it ideal for daily lab routine
- Easy and quick fixation of denudation pipettes with different sizes

Ordering Information

Type	Ref Nr.	Material	Qty/Unit	
Handset	V-DH02	Metal	1	





V-SPERM MORPHO SLIDE is intended for staining of human sperm and qualitative assessment of sperm morphology for in vitro diagnostic purposes in IVF/ART.

Details

- Stains: cresyl violet acetate and new methyl blue
- Time and cost effecient
- Ready to use
- 1 x 10 pre-stained slides and 10 coverslips

Ordering Information

Туре	Ref Nr.	Stains	Qty/Unit
Pre-stained Slides	V-SMS-01	Cresyl violet acetate & new methylene blue	10 (Slides and Cover slides)



V-GLASSPET pipettes are designed for easy liquid handling.

Details

- Endotoxin tested (< 0.05 EU/device)
- MEA 1-Cell to expanded blastocyst in 96h: ≥ 80%
- Steam sterilized, available plugged and non-plugged
- Sizes: 150 mm and 230 mm
- 2 mL volume
- 5 pcs/pouch, 5 pouches/box, 10 boxes/unit → 250 pcs/unit

Ordering Information

Type	Ref Nr.	Dimension (mm)	Classification	Qty/Unit
Pasteur Pipette	V-GPT-150	150 (Length)	Plugged	250
Pasteur Pipette	V-GPT-150NP	150 (Length)	Non-Plugged	250
Pasteur Pipette	V-GPT-230	230 (Length)	Plugged	250
Pasteur Pipette	V-GPT-230NP	230 (Length)	Non-Plugged	250

MARKER PEN



V-MARKER PEN with its ultra fine tip are ideal for the marking of dishes, cryovials and straws.

Details

- Endotoxin tested (< 0.05 EU/device)
- MEA 1-Cell to expanded blastocyst in 96h: ≥ 80%
- Non-toxic and xylene free
- Available in black, red and blue
- 1 x 4 pens in your choice of colors

Ordering Information

Туре	Ref Nr.	Classification	Qty/Unit
Marker Pen	V-PEN-DC	Red, Blue & Black	4
Marker Pen	V-PEN-BK	Black	4
Marker Pen	V-PEN-BL	Blue	4
Marker Pen	V-PEN-R	Red	4

PENTOXIFYLLINE



V-PENTOXIFYLLINE is a HEPES buffered HTF medium enriched with pentoxifillyine. It increases the short-term motility of sperm cells enabling the differentiation between dead and immotile sperm originating from TESE or necrozoospermic samples.

Details

- Endotoxin tested (< 0.5 EU/mL)
- Store at 2-8 °C
- 5 x 1mL Microvial
- Human Sperm Survival Assay tested ≥ 80%

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
For sperm activation	V-PENTO	1	PPCO	5×1

for use in reasearch procedures only

YONOPHORE



 $V-IONOPHORE\ induces\ Ca-dependent\ Oocyte\ activation\ improving\ fertilization\ outcomes\ especially\ with\ prior\ fertilization\ failures\ using\ ICSI\ procedures.$

Details

- Endotoxin tested (< 0.25 EU/mL) 5x 1mL
- MEA 1-Cell to expanded blastocyst in 96h: ≥ 80%
- Store at 2-8 °C
- 5 x 1mL Microvial

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
For Oocyte activation	V-IONO	1	PPCO	5 x 1

for use in reasearch procedures only

Info

Please check the product availability in your region while ordering.





V-PURE is sterile water.

Details

- Endotoxin tested
- 1 x 2000 mL PET bottle

Ordering Information

Type	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
Sterile Water	V-PW-2000	2000	PET	1

DENUPET ERGOGRIP

 $V-DENUPET\ ErgoGrip\ is\ an\ accessory\ for\ V-DENUPET\ HANDLE\ DH02.\ It's\ made\ of\ medical\ grade\ ABS.\ The\ ErgoGrip\ gives\ comfort\ while\ handling\ the\ DENUPET\ HANDLE.\ The\ ergonomic\ design\ reduces\ the\ stress\ on\ the\ fingers\ during\ the\ usage.$

Details

• Made from easy-to-clean ABS material



Туре	Ref Nr.	Qty/Unit
ErgoGrip	V-DH02-EG	1

GLASSPET BULBS



Silicone pipette bulb to use with any pipettes for easy liquid handling.

Details

- Made from biocompatible silicone
- Available in 5 colour variants
- Plastic ziplock bags with labelling with 25 pcs inside (5 pieces of each colour)

Ordering Information

Туре	Ref Nr.	Classification	Qty/Unit
Pipetting Bulb	V-RUT-MIX	5 colour variants (Yellow, Pink, Red, Turquoise, Blue)	25 (5 pcs / each colour)

HANDPET



Plastic pipet for liquid handling

Details

- Made of polycarbonate and PP/TPE
- Individual soft blister packed
- Sterilized by radiation
- Compatible with luer lock syringes

Ordering Information

Туре	Ref Nr.	Dimentions	Qty/Unit	
Plastic hand pipet	V-HPT-90	90	250	

6 WELLDISH



V-6 WELLDISH for laboratory use

Details

- 6 well dish with flat bottom
- 6 well dish with flat bottom

 Clear, untreated virgin polystyrene with excellent optical properties
 Size: 64 x 77 x 10 mm

 Volume of each well: 0.25 mL

 Diameter of each well: 15 mm

 200 pcs/box

 Single soft blister packed

Ordering Information

Туре	Ref Nr.	Dimentions	Qty/Unit	
6-well-dish	V-6WD	64 x 77 x 10 mm	200	

V-CARE







 $V\text{-}CARE\ STERI\ CONTROL\ enables\ an\ easy\ control\ of\ the\ aseptic\ working\ environment\ and\ medium\ handling. The\ product\ enables\ an\ unmatched\ level\ of\ control\ over\ the\ sterile\ liquid\ handling\ processes\ used\ in\ ART.\ V\text{-}CARE\ SteriControl\ ensures\ the\ traceability\ of\ sterile\ working\ conditions.$

Details

- Tryptic Soy Broth (TSB) medium
- Growth promotion tested
- Shelf life 12 months from date of production

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
STERI CONTROL	V-CARE-SC20P	20	PETG	1
STERI CONTROL	V-CARE-SC20G	20	GLASS	1





V-CARE INCU CHECK enables high-volume pH measurement, designed for frequent pH monitoring, directly linked to incubator CO_2 levels.

Details

- CO₂ dependent buffered solution
- Sterility tested
- Shelf life 12 months from date of production

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
INCU CHECK	V-CARE-IC20G	20	GLASS	1





 $V\text{-}\mathsf{CARE}\ \mathsf{CLEAN}\ \mathsf{CHECK}\ \mathsf{helps}\ \mathsf{evaluate}\ \mathsf{the}\ \mathsf{thoroughness}\ \mathsf{of}\ \mathsf{cleaning}\ \mathsf{on}\ \mathsf{embryology/andrology}\ \mathsf{lab}\ \mathsf{and}\ \mathsf{other}\ \mathsf{critical}\ \mathsf{surfaces}.$

Details

- Fast drying substance
- Human Sperm Survival Assay tested ≥ 80% Control Motility
- Shelf life 12 months from date of production
- 1 x 0,8 mL Microvial with 3 Cotton Swabs included

Ordering Information

Туре	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
CLEAN CHECK	V-CARE-CCP	0.8	Cryovial	1



Our main offices



JENA

Göschwitzer Str. 22 07745 Jena Germany Phone: +49 3641 539 19 76



LANGENFELD

Raiffeisenstr. 15a 40764 Langenfeld Germany Phone: +49 2173 200 41 30



Scan me

Connect with VITROMEDonline







